

Sex hormone-binding globulin

Catalog No: #AP89477



Package Size: #AP89477-1 50ug #AP89477-2 100ug #AP89477-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sex hormone-binding globulin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	SHBG, Sex steroid-binding protein, SBP, Testis-specific androgen-binding protein, ABP
Accession No.	Uniprot ID:P97497
Uniprot	P97497
GeneID	20415;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Shbg

Note: This product is for in vitro research use only