

SPARC-like protein 1

Catalog No: #AP89483



Package Size: #AP89483-1 50ug #AP89483-2 100ug #AP89483-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	SPARC-like protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	High endothelial venule protein,Hevin,MAST 9
Accession No.	Uniprot ID:Q14515
Uniprot	Q14515
GeneID	8404;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SPARCL1

Note: This product is for in vitro research use only