

Stanniocalcin-1

Catalog No: #AP89485



Package Size: #AP89485-1 50ug #AP89485-2 100ug #AP89485-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Stanniocalcin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	STC-1,STC
Accession No.	Uniprot ID:Q9N0T1
Uniprot	Q9N0T1
GeneID	338078;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:STC1

Note: This product is for in vitro research use only