

# Alpha-synuclein

Catalog No: #AP89486



Package Size: #AP89486-1 50ug #AP89486-2 100ug #AP89486-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Alpha-synuclein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Non-A beta component of AD amyloid,Non-A4 component of amyloid precursor,NACP,Syn
Accession No.	Uniprot ID:O55042
Uniprot	O55042
GeneID	20617;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Snca

Note: This product is for in vitro research use only