

Transcription factor p65

Catalog No: #AP89489



Package Size: #AP89489-1 50ug #AP89489-2 100ug #AP89489-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Transcription factor p65
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Nuclear factor NF-kappa-B p65 subunit
Accession No.	Uniprot ID:P98152
Uniprot	P98152
GeneID	396027;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RELA

Note: This product is for in vitro research use only