

## RXFP1 Antibody

Catalog No: #44621

Package Size: #44621-1 50ul #44621-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

|                       |   |
|-----------------------|---|
| Product Name          | RXFP1 Antibody  |
| Host Species          | Rabbit  |
| Clonality             | Polyclonal  |
| Purification          | Immunogen affinity purified   |
| Applications          | WB  |
| Species Reactivity    | Hu Ms Rt  |
| Specificity           | The antibody was purified by immunogen affinity chromatography.   |
| Immunogen Type        | Peptide   |
| Immunogen Description | Recombinant fusion protein of human LGR7  |
| Target Name           | RXFP1   |
| Other Names           | RXFP1, LGR7, LGR7.1, LGR7.10, LGR7.2, RXFPR1, Relaxin receptor 1  |
| Accession No.         | Swiss-Prot#: Q9HBX9NCBI Gene ID:59350   |
| Uniprot               | Q9HBX9  |
| GeneID                | 59350;  |
| Calculated MW         | Predicted band size: 21; 75; 81; 83; 86 kD; Observed band size: 87 kD                                     |
| Concentration         | 1.0mg/mL  |
| Formulation           | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide. |
| Storage               | Store at -20°C  |

## Application Details

WB 1:500-2000

## Images



Western blot analysis of LGR7 expression in HeLa (A), mouse brain (B), rat testis (C) whole cell lysates.

---

Note: This product is for in vitro research use only