

Transforming growth factor beta receptor type 3

Catalog No: #AP89490



Package Size: #AP89490-1 50ug #AP89490-2 100ug #AP89490-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Transforming growth factor beta receptor type 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	TGF-beta receptor type 3, Betaglycan, Transforming growth factor beta receptor III, TGF-beta receptor type III
Accession No.	Uniprot ID:Q03167
Uniprot	Q03167
GeneID	7049;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TGFBR3

Note: This product is for in vitro research use only