

Thromboxane-A synthase

Catalog No: #AP89491



Package Size: #AP89491-1 50ug #AP89491-2 100ug #AP89491-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Thromboxane-A synthase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	TS,Cytochrome P450 5A1,Cyp5,Cyp5a1
Accession No.	Uniprot ID:P36423
Uniprot	P36423
GeneID	21391;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tbxas1

Note: This product is for in vitro research use only