

## Thimet oligopeptidase

Catalog No: #AP89492



Package Size: #AP89492-1 50ug #AP89492-2 100ug #AP89492-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Thimet oligopeptidase  |
| Brief Description  | Recombinant Protein  |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE   |
| Species Reactivity | Rt   |
| Other Names        | Endo-oligopeptidase A,Endopeptidase 24.15,PZ-peptidase,Soluble metallo-endopeptidase |
| Accession No.      | Uniprot ID:P24155  |
| Uniprot            | P24155   |
| GeneID             | 64517;   |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)                               |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Thop1

Note: This product is for in vitro research use only