

Metalloproteinase inhibitor 3

Catalog No: #AP89494



Package Size: #AP89494-1 50ug #AP89494-2 100ug #AP89494-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Metalloproteinase inhibitor 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Tissue inhibitor of metalloproteinases 3, TIMP-3, Sun, Timp-3
Accession No.	Uniprot ID:P39876
Uniprot	P39876
GeneID	21859;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Timp3

Note: This product is for in vitro research use only