

Tumor necrosis factor receptor superfamily member 6

Catalog No: #AP89497



Package Size: #AP89497-1 50ug #AP89497-2 100ug #AP89497-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor receptor superfamily member 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Apo-1 antigen, Apoptosis-mediating surface antigen FAS, FASLG receptor, CD95, Apt1, Tnfrsf6
Accession No.	Uniprot ID:Q63199
Uniprot	Q63199
GeneID	246097;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Fas

Note: This product is for in vitro research use only