

Tryptase delta

Catalog No: #AP89499



Package Size: #AP89499-1 50ug #AP89499-2 100ug #AP89499-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tryptase delta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Delta-tryptase,HmMCP-3-like tryptase III,Mast cell mMCP-7-like,Tryptase-3
Accession No.	Uniprot ID:Q9BZJ3
Uniprot	Q9BZJ3
GeneID	23430;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TPSD1

Note: This product is for in vitro research use only