

Thymidine phosphorylase

Catalog No: #AP89501



Package Size: #AP89501-1 50ug #AP89501-2 100ug #AP89501-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Thymidine phosphorylase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	TP,TdRPase,Ecgf1
Accession No.	Uniprot ID:Q5FVR2
Uniprot	Q5FVR2
GeneID	315219;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tymp

Note: This product is for in vitro research use only