

WNT1-inducible-signaling pathway protein 3

Catalog No: #AP89505



Package Size: #AP89505-1 50ug #AP89505-2 100ug #AP89505-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	WNT1-inducible-signaling pathway protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	WISP-3,CCN family member 6,CCN6,UNQ462/PRO790/PRO956
Accession No.	Uniprot ID:O95389
Uniprot	O95389
GeneID	8838;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:WISP3

Note: This product is for in vitro research use only