

Cytidine deaminase

Catalog No: #AP89507



Package Size: #AP89507-1 50ug #AP89507-2 100ug #AP89507-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cytidine deaminase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cytidine aminohydrolase,Cdd
Accession No.	Uniprot ID:P56389
Uniprot	P56389
GeneID	72269;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cda

Note: This product is for in vitro research use only