

Glutamyl aminopeptidase

Catalog No: #AP89509



Package Size: #AP89509-1 50ug #AP89509-2 100ug #AP89509-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glutamyl aminopeptidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	EAP,Aminopeptidase A,AP-A,CD249
Accession No.	Uniprot ID:Q32LQ0
Uniprot	Q32LQ0
GeneID	504350;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ENPEP

Note: This product is for in vitro research use only