

C-X-C motif chemokine 6

Catalog No: #AP89510

Package Size: #AP89510-1 50ug #AP89510-2 100ug #AP89510-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-X-C motif chemokine 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Granulocyte chemotactic protein 2,GCP-2,Small-inducible cytokine B6,GCP2,SCYB6
Accession No.	Uniprot ID:P80221
Uniprot	P80221
GeneID	281735;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CXCL6

Note: This product is for in vitro research use only