

Bovine ALCAM/CD166 antigen ELISA Kit

Catalog No: #EK19155

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine ALCAM/CD166 antigen ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine CD166 antigen
Target Name	CD166 antigen
Other Names	Activated leukocyte cell adhesion molecule, CD166
Accession No.	Q9BH13
Uniprot	Q9BH13
GeneID	281614;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.078ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: ALCAM

Note: This product is for in vitro research use only