

## Dog ALCAM/CD166 antigen ELISA Kit

Catalog No: #EK19157

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Dog ALCAM/CD166 antigen ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Dog
Specificity	Natural and recombinant Dog CD166 antigen
Target Name	CD166 antigen
Other Names	Activated leukocyte cell adhesion molecule, SB-10 antigen, CD166
Accession No.	O46634
Uniprot	O46634

## Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.078ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: ALCAM

Note: This product is for in vitro research use only