

Rat Alcam/CD166 antigen ELISA Kit

Catalog No: #EK19158

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Alcam/CD166 antigen ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat CD166 antigen
Target Name	CD166 antigen
Other Names	Activated leukocyte cell adhesion molecule, HB2, KG-CAM, Protein MEMD, SB-10 antigen, CD166
Accession No.	O35112
Uniprot	O35112
GeneID	79559;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.078ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Alcam

Note: This product is for in vitro research use only