

## Rat Areg/Amphiregulin ELISA Kit

Catalog No: #EK19160

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Areg/Amphiregulin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Amphiregulin
Target Name	Amphiregulin
Other Names	AR,Schwannoma-derived growth factor,SDGF,SdGF
Accession No.	P24338
Uniprot	P24338
GeneID	29183;

## Application Details

Assay Type: Sandwich Detect Range:15.6-1000pg/mL Sensitivity:7.83pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Areg

Note: This product is for in vitro research use only