

Chicken ANG/Angiogenin ELISA Kit

Catalog No: #EK19161

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Chicken ANG/Angiogenin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Chicken
Specificity	Natural and recombinant Chicken Angiogenin
Target Name	Angiogenin
Other Names	Ribonuclease A
Accession No.	P27043
Uniprot	P27043
GeneID	396194;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: ANG

Note: This product is for in vitro research use only