

Bovine ANGPT1/Angiopoietin-1 ELISA Kit

Catalog No: #EK19163

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine ANGPT1/Angiopoietin-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Angiopoietin-1
Target Name	Angiopoietin-1
Other Names	ANG-1,ANG1
Accession No.	O18920
Uniprot	O18920
GeneID	282140;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.15ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:ANGPT1

Note: This product is for in vitro research use only