

Dog ANGPT1/Angiopoietin-1 ELISA Kit

Catalog No: #EK19164

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Dog ANGPT1/Angiopoietin-1 ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Dog |
| Specificity | Natural and recombinant Dog Angiopoietin-1 |
| Target Name | Angiopoietin-1 |
| Other Names | ANG-1 |
| Accession No. | Q60FC1 |
| Uniprot | Q60FC1 |
| GeneID | 403656; |

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:ANGPT1

Note: This product is for in vitro research use only