

Pig ANGPT1/Angiopoietin-1 ELISA Kit

Catalog No: #EK19165

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pig ANGPT1/Angiopoietin-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig Angiopoietin-1
Target Name	Angiopoietin-1
Other Names	ANG-1
Accession No.	Q9BDY8
Uniprot	Q9BDY8
GeneID	397009;

Application Details

Assay Type: Sandwich Detect Range: 3.12-200ng/mL Sensitivity: 1.22ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: ANGPT1

Note: This product is for in vitro research use only