

## Bovine ANGPT2/Angiopoietin-2 ELISA Kit

Catalog No: #EK19166

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bovine ANGPT2/Angiopoietin-2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Angiopoietin-2
Target Name	Angiopoietin-2
Other Names	ANG-2,ANG2
Accession No.	O77802
Uniprot	O77802
GeneID	282141;

## Application Details

Assay Type: Sandwich Detect Range:78-5000pg/ml Sensitivity:38pg/ml Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:ANGPT2

Note: This product is for in vitro research use only