

# Bovine CTGF/Connective tissue growth factor ELISA Kit



Catalog No: #EK19167

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bovine CTGF/Connective tissue growth factor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Connective tissue growth factor
Target Name	Connective tissue growth factor
Other Names	CCN family member 2,CCN2
Accession No.	O18739
Uniprot	O18739
GeneID	281103;

## Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.15ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CTGF

Note: This product is for in vitro research use only