

Pig CTGF/Connective tissue growth factor ELISA Kit

Catalog No: #EK19168

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pig CTGF/Connective tissue growth factor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig Connective tissue growth factor
Target Name	Connective tissue growth factor
Other Names	CCN family member 2,CCN2
Accession No.	O19113
Uniprot	O19113

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.15ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CTGF

Note: This product is for in vitro research use only