

Bovine BDNF/Brain-derived neurotrophic factor ELISA Kit



Catalog No: #EK19169

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine BDNF/Brain-derived neurotrophic factor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Brain-derived neurotrophic factor
Target Name	Brain-derived neurotrophic factor
Other Names	BDNF
Accession No.	Q95106
Uniprot	Q95106
GeneID	617701;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:BDNF

Note: This product is for in vitro research use only