

Bovine BTC/Probetacellulin ELISA Kit

Catalog No: #EK19170



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bovine BTC/Probetacellulin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Probetacellulin
Target Name	Probetacellulin
Accession No.	Q9TTC5
Uniprot	Q9TTC5

Application Details

Assay Type: Sandwich Detect Range: 62.5-4000pg/mL Sensitivity: 31.2pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: BTC

Note: This product is for in vitro research use only