

Mouse Btc/Probetacellulin ELISA Kit

Catalog No: #EK19171

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Btc/Probetacellulin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Probetacellulin
Target Name	Probetacellulin
Other Names	Bcn
Accession No.	Q05928
Uniprot	Q05928
GeneID	12223;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.081ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Btc

Note: This product is for in vitro research use only