

Chicken BMP2/Bone morphogenetic protein 2 ELISA Kit



Catalog No: #EK19172

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Chicken BMP2/Bone morphogenetic protein 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Chicken
Specificity	Natural and recombinant Chicken Bone morphogenetic protein 2
Target Name	Bone morphogenetic protein 2
Other Names	BMP-2,BMP-2
Accession No.	Q90751
Uniprot	Q90751

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:BMP2

Note: This product is for in vitro research use only