

Rat Bmp2/Bone morphogenetic protein 2 ELISA Kit

Catalog No: #EK19173

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Bmp2/Bone morphogenetic protein 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Bone morphogenetic protein 2
Target Name	Bone morphogenetic protein 2
Other Names	BMP-2, Bone morphogenetic protein 2A, BMP-2A, Bmp-2
Accession No.	P49001
Uniprot	P49001
GeneID	29373;

Application Details

Assay Type: Sandwich Detect Range: 62.5-4000pg/mL Sensitivity: 31.25pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Bmp2

Note: This product is for in vitro research use only