

Bovine BMP4/Bone morphogenetic protein 4 ELISA Kit



Catalog No: #EK19174

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine BMP4/Bone morphogenetic protein 4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Bone morphogenetic protein 4
Target Name	Bone morphogenetic protein 4
Other Names	BMP-4
Accession No.	Q2KJH1
Uniprot	Q2KJH1
GeneID	407216;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 37pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: BMP4

Note: This product is for in vitro research use only