

Chicken BMP4/Bone morphogenetic protein 4 ELISA Kit



Catalog No: #EK19175

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Chicken BMP4/Bone morphogenetic protein 4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Chicken
Specificity	Natural and recombinant Chicken Bone morphogenetic protein 4
Target Name	Bone morphogenetic protein 4
Other Names	BMP-4,BMP-4
Accession No.	Q90752
Uniprot	Q90752

Application Details

Assay Type: SandwichDetect Range:15.6-1000pg/mL Sensitivity:5.4pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:BMP4

Note: This product is for in vitro research use only