

# Mouse Bmpr1a/Bone morphogenetic protein receptor type-1A ELISA Kit

Catalog No: #EK19177

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Bmpr1a/Bone morphogenetic protein receptor type-1A ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Bone morphogenetic protein receptor type-1A
Target Name	Bone morphogenetic protein receptor type-1A
Other Names	BMP type-1A receptor,Activin receptor-like kinase 3,ALK-3,BMP-2/BMP-4 receptor,Serine/threonine-protein kinase receptor R5,SKR5,CD292,Acvrk3,Bmpr
Accession No.	P36895
Uniprot	P36895
GeneID	12166;

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.158ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Bmpr1a

Note: This product is for in vitro research use only