

Human BMPR2/Bone morphogenetic protein receptor type-2 ELISA Kit



Catalog No: #EK19178

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human BMPR2/Bone morphogenetic protein receptor type-2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Bone morphogenetic protein receptor type-2
Target Name	Bone morphogenetic protein receptor type-2
Other Names	BMP type-2 receptor,Bone morphogenetic protein receptor type II,BMP type II receptor,PPH1
Accession No.	Q13873
Uniprot	Q13873
GeneID	659;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.082ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:BMPR2

Note: This product is for in vitro research use only