

Mouse Bmpr2/Bone morphogenetic protein receptor type-2 ELISA Kit

Catalog No: #EK19179

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Bmpr2/Bone morphogenetic protein receptor type-2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Bone morphogenetic protein receptor type-2
Target Name	Bone morphogenetic protein receptor type-2
Other Names	BMP type-2 receptor, BRK-3, Bone morphogenetic protein receptor type II, BMP type II receptor
Accession No.	O35607
Uniprot	O35607
GeneID	12168;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.083ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Bmpr2

Note: This product is for in vitro research use only