

## Bovine CDH1/Cadherin-1 ELISA Kit

Catalog No: #EK19180

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bovine CDH1/Cadherin-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Cadherin-1
Target Name	Cadherin-1
Other Names	Epithelial cadherin,E-cadherin,CD324
Accession No.	Q6R8F2
Uniprot	Q6R8F2
GeneID	282637;

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:CDH1

Note: This product is for in vitro research use only