

Chicken CDH1/Cadherin-1 ELISA Kit

Catalog No: #EK19181

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Chicken CDH1/Cadherin-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Chicken
Specificity	Natural and recombinant Chicken Cadherin-1
Target Name	Cadherin-1
Other Names	Epithelial cadherin,E-cadherin,Liver cell adhesion molecule,L-CAM
Accession No.	P08641
Uniprot	P08641
GeneID	415860;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CDH1

Note: This product is for in vitro research use only