

Dog CDH1/Cadherin-1 ELISA Kit

Catalog No: #EK19182

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dog CDH1/Cadherin-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Dog
Specificity	Natural and recombinant Dog Cadherin-1
Target Name	Cadherin-1
Other Names	Epithelial cadherin,E-cadherin,Uvomorulin,CD324
Accession No.	F1PAA9
Uniprot	F1PAA9
GeneID	442858;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/ml Sensitivity:0.052ng/ml Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CDH1

Note: This product is for in vitro research use only