

## Mouse Cdh1/Cadherin-1 ELISA Kit

Catalog No: #EK19183

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Cdh1/Cadherin-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cadherin-1
Target Name	Cadherin-1
Other Names	ARC-1, Epithelial cadherin, E-cadherin, Uvomorulin, CD324
Accession No.	P09803
Uniprot	P09803
GeneID	12550;

## Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39.1pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Cdh1

Note: This product is for in vitro research use only