

Rat Cdh1/Cadherin-1 ELISA Kit

Catalog No: #EK19184

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Cdh1/Cadherin-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Cadherin-1
Target Name	Cadherin-1
Other Names	Epithelial cadherin,E-cadherin,Uvomorulin,CD324
Accession No.	Q9R0T4
Uniprot	Q9R0T4
GeneID	83502;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.084ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Cdh1

Note: This product is for in vitro research use only