

Mouse Cntf/Ciliary neurotrophic factor ELISA Kit

Catalog No: #EK19186

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cntf/Ciliary neurotrophic factor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Ciliary neurotrophic factor
Target Name	Ciliary neurotrophic factor
Other Names	CNTF
Accession No.	P51642
Uniprot	P51642
GeneID	12803;

Application Details

Assay Type: Sandwich Detect Range: 39-2500pg/mL Sensitivity: 19.5pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cntf

Note: This product is for in vitro research use only