

Pig CNTF/Ciliary neurotrophic factor ELISA Kit

Catalog No: #EK19187

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pig CNTF/Ciliary neurotrophic factor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig Ciliary neurotrophic factor
Target Name	Ciliary neurotrophic factor
Other Names	CNTF
Accession No.	O02732
Uniprot	O02732

Application Details

Assay Type: Sandwich Detect Range: 31.2-2000pg/mL Sensitivity: 15.6pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CNTF

Note: This product is for in vitro research use only