

## Bovine CORO1A/Coronin-1A ELISA Kit

Catalog No: #EK19188

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bovine CORO1A/Coronin-1A ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Coronin-1A
Target Name	Coronin-1A
Other Names	Coronin-like protein A,Clipin-A,Coronin-like protein p57,Tryptophan aspartate-containing coat protein,TACO,CORO1
Accession No.	Q92176
Uniprot	Q92176
GeneID	282196;

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:CORO1A

Note: This product is for in vitro research use only