

Mouse Coro1a/Coronin-1A ELISA Kit

Catalog No: #EK19189

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Coro1a/Coronin-1A ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Coronin-1A
Target Name	Coronin-1A
Other Names	Coronin-like protein A,Clipin-A,Coronin-like protein p57,Tryptophan aspartate-containing coat protein,TACO,Coro1
Accession No.	O89053
Uniprot	O89053
GeneID	12721;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.086ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Coro1a

Note: This product is for in vitro research use only