

Rat Coro1a/Coronin-1A ELISA Kit

Catalog No: #EK19190

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Coro1a/Coronin-1A ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Coronin-1A
Target Name	Coronin-1A
Other Names	Coronin-like protein A,Clipin-A,Tryptophan aspartate-containing coat protein,TACO,Coro1
Accession No.	Q91ZN1
Uniprot	Q91ZN1
GeneID	155151;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.087ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Coro1a

Note: This product is for in vitro research use only