

Bovine CSF3/Granulocyte colony-stimulating factor ELISA Kit



Catalog No: #EK19191

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine CSF3/Granulocyte colony-stimulating factor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Granulocyte colony-stimulating factor
Target Name	Granulocyte colony-stimulating factor
Other Names	G-CSF,GCSF
Accession No.	P35833
Uniprot	P35833
GeneID	281096;

Application Details

Assay Type: SandwichDetect Range:78-5000pg/mL Sensitivity:39pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:CSF3

Note: This product is for in vitro research use only