

Dog CSF3/Granulocyte colony-stimulating factor ELISA Kit



Catalog No: #EK19192

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dog CSF3/Granulocyte colony-stimulating factor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Dog
Specificity	Natural and recombinant Dog Granulocyte colony-stimulating factor
Target Name	Granulocyte colony-stimulating factor
Other Names	G-CSF
Accession No.	P35834
Uniprot	P35834

Application Details

Assay Type: Sandwich Detect Range: 15.6-1000pg/mL Sensitivity: 4.6pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CSF3

Note: This product is for in vitro research use only